T.I.T.L. MARKET SAMPLE NO. 4

TEST 1 AND 2 RESULTS

GROUPED BY SMOKING MACHINE

		<u>"Tar"</u>	mg/cigt Te	est 1	<u>"Tar</u>	<u>"Tar" mg/cigt Tes</u>		
BRAND		Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	
Alpine	X s CV	19.4 ± 0.3 0.5 2.7	19.4±0.8 1.2 6.2	19.4±0.4 0.9 4.6	19.4±0.4 0.6 3.1	19.0±0.5 0.8 4.5	19.2±0.3 0.7 3.9	
Belair	x s cv	$ \begin{array}{c} 18.6 \pm 0.6 \\ 1.0 \\ 5.2 \end{array} $	18.9 ± 0.6 0.9 4.8	18.8±0.4 0.9 5.0	18.0±0.3 0.5 2.8	17.4±0.4 0.7 3.9	17.7 ± 0.3 0.6 3.6	
Benson & Hedges reg. size	X S CV	$19.6 \pm 0.7 \\ 1.2 \\ 6.0$	20.4 ± 0.7 1.1 5.3	20.0 ± 0.5 1.2 5.8	15.9 ± 0.7 1.1 6.7	$ \begin{array}{r} 15.4 \pm 0.3 \\ 0.4 \\ 2.5 \end{array} $	15.7±0.4 0.8 5.2	
Benson & Hedges king size	X S CV	20.2 ± 0.3 0.4 2.2	20.0±0.5 0.8 4.0	$\begin{array}{c} 20.1 \pm 0.3 \\ 0.6 \\ 3.2 \end{array}$	18.2 ± 0.5 0.8 4.4	18.0±0.6 0.9 5.0	18.1±0.4 0.8 4.6	
Benson & Hedges 100 mm. filter	X s CV	21.4±0.4 0.7 3.1	20.9 ± 0.5 0.8 3.7	21.2±0.3 0.7 3.5	21.6±0.5 0.8 3.5	21.1±0.5 0.8 3.9	21.4±0.4 0.8 3.8	
Benson & Hedges 100 mm. menthol	x s CV	22.1 ± 0.4 0.7 3.3	22.1±0.6 0.9 4.2	22.1±0.4 0.8 3.6	22.0 ± 0.4 0.7 3.2	22.1±0.4 0.6 2.8	22.0±0.3 0.6 2.9	
Bull Durham	CA a X	29.2±0.5 0.8 2.8	29.2 ± 0.5 0.8 2.7	29.2±0.3 0.8 2.6	26.9 ± 0.5 0.8 3.1	26.3±0.6 0.9 3.5	26.6±0.4 0.9 3.4	

		"Tar	'mg/cigt Te	est 1	''Tar	" mg/cigt Te	est 2
BRAND		Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined
Camel	X	26.4±0.8	25.7 ± 0.8	26.1±0.6	21.4 ± 0.5	21.0 ± 0.7	21.2± 0.4
reg. size,	s	1.2	1.3	1.3	0.8	1.2	1.0
non-filter	CV	4.5	5.0	4.8	4.0	5.6	4.8
Camel	X	20.5 ± 0.5	20.8 ± 0.6	20.7 ± 0.4	19.4 ± 0.4	18.9 ±0.3	19.2 ±0.3
king size	s	0.8	1.0	0.9	0.7	0.5	0.6
filter	CV	3.7	4.7	4.2	3.4	2.6	3.3
Carlton	X S CV	4.0±0.6 1.0 25.0	4.1±0.6 1.0 24.4	4.0±0.4 1.0 24.1	4.5 ± 0.4 0.6 13.6	4.3 ±0.5 0.7 17.0	4.4 ±0.3 0.7 15.2
Cascade	x s cv	7.8 ± 0.5 0.8 10.2	7.1 ± 0.5 0.8 12.0	7.5±0.4 0.9 11.7	6.6 ± 0.2 0.4 5.9	6.2±0.4 0.6 10.1	6.4±0.2 0.5 8.5
Century	X	22.7 ± 0.6	22.7 ± 0.5	22.7±0.4	22.8 ± 0.3	22.6 ±0.5	22.7±0.3
	S	1.0	0.8	0.9	0.4	0.8	0.6
	CV	4.3	3.8	3.9	2.0	3.6	2.9
Chesterfield	X	24.8 ± 0.7	24.9±0.4	24.9±0.4	20.6 ± 0.6 1.0 4.7	19.7±0.2	20.1±0.4
reg. size	s	1.1	0.6	0.9		0.4	0.8
non-filter	CV	4.5	2.6	3.6		1.8	4.2
Chesterfield	X	29.3±0.7	29.4±0.9	$ \begin{array}{c} 29.4 \pm 0.5 \\ 1.2 \\ 4.2 \end{array} $	25.1±0.6	24.2±0.5	24.7±0.4
king size	S	1.0	1.4		0.9	0.7	0.9
non-filter	CV	3.6	4.8		3.6	3.0	3.7
Chesterfield king size filter	X s CV	22.2 ± 0.8 1.2 5.4	22.0 ± 0.5 0.8 3.7	22.1 ± 0.4 1.0 4.6	20.8±0.3 0.5 2.6	20.1±0.4 0.6 ·3.4	20.4±0.3 0.7 3.4

υ
_
_
ω
0
0
0
9
Ω
N
\prec
\rightharpoonup
7

•		"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	
Chesterfield king size, filter menthol	X s CV	20.2 ± 1.0 1.6 7.8	19.6 ± 0.8 1.2 6.4	19.9± 0.6 1.4 7.1	18.4 ± 0.4 0.7 3.6	18.2±0.4 0.6 3.1	18.3±0.3 0.6 3.3	
Chesterfield 101 mm.	X S CV	20.2 ± 0.6 1.0 5.0	19.6 ± 0.6 1.0 5.3	20.0±0.5 1.0 5.2	$\begin{array}{c} 20.4 \pm 0.4 \\ 0.7 \\ 3.2 \end{array}$	20.0±0.8 1.3 6.4	20.2 ± 0.4 1.0 5.1	
Colony 100 mm. filter	X s CV	28.2 ± 0.7 1.1 3.9	28.0 ± 1.0 1.6 5.8	28.0±0.6 1.4 4.8	26.2 ± 0.4 0.6 2.4	26.3 ± 0.4 0.6 2.5	26.2±0.3 .0.6 2.4	
Colony 100 mm. filter menthol	X s CV	25.7 ± 0.6 0.9 3.6	25.9±0.8 1.2 4.7	25.8±0.5 1.0 4.1	$ \begin{array}{c} 24.1 \pm 0.4 \\ 0.7 \\ 2.9 \end{array} $	23.6 ± 0.5 0.8 3.6	23.9±0.4 0.8 3.4	
Domino king size, non-filter	X s CV	25.4 ± 1.2 2.0 7.7	25.9 ± 0.5 0.8 3.3	25.7±0.7 1.5 5.8	22.5 ± 0.8 1.2 5.4	22.0±0.5 0.8 3.8	22.2±0.5 1.0 4.7	
Domino king size, filter	X s CV	20.3±0.6 1.0 4.9	21.2 ± 0.5 0.8 3.5	20.8±0.4 1.0 4.6	17.5±0.4 0.6 3.5	17.2±0.6 0.9 5.4	17.4±0.4 0.8 4.6	
Domino king size, filter menthol	X s CV	13.6±1.0 1.6 11.6	15.0 ± 0.8 1.3 8.5	14.3 ± 0.7 1.6 11.0	11.4 ± 0.5 0.7 6.6	$11.1 \pm 0.7 \\ 1.2 \\ 10.5$	11.3±0.4 1.0 8.6	
Duke	X s CV	9.7 ± 0.2 0.4 4.2	9.4±0.6 0.9 9.3	9.6 ±0.3 0.7 7.1	9.9±0.6 0.9 9.4	9.9 ± 0.4 0.6 6.2	9.9 ±0.3 0.8 7.8	

		"Tar"	mg/cigt Test	<u> </u>	"Tar"	mg/cigt Tes	eigt Test 2	
BRAND	,	Smoking <u>Machine A</u>	Smoking <u>Machine</u> B	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine</u> B	A & B Combined	
Du Maurier	x s CV	18.3±0.6 1.0 5.5	18.2±0.3 0.5 2.9	18.2 ± 0.4 0.8 4.3	17.3 ± 0.5 0.7 4.2	17.2 ± 0.2 0.4 2.3	17.2 ± 0.2 0.6 3.3	
English Ovals	X	27.3±0.8	26.6±1.3	27.0 ±0.7	22.6 ± 0.7 1.1 5.0	21.8±0.3	22.2±0.4	
reg. size,	s	1.3	2.0	1.7		0.4	0.9	
non-filter	CV	4.6	7.5	6.2		1.9	4.3	
English Ovals king size, non-filter	x	35.0±1.3	35.3±0.8	35.1 ±0.8	29.4 ± 0.7	30.4±0.8	29.9±0.6	
	s	2.0	1.3	1.7	1.1	1.3	1.3	
	CV	5.9	3.7	4.8	3.7	4.3	4.3	
Fatima	x s CV	33.3± 0.8 1.2 3.6	33.1±0.7 1.1 3.4	33.2±0.5 1.1 3.4	28.3±1.1 1.8 6.3	27.4 ± 0.7 1.1 4.2	27.8±0.7 1.5 5.5	
Frappe	X	9.0±0.9	9.0 ±0.6	9.0±0.5	9.0±1.1	7.8±0.5	8.4±0.6	
	s	1.4	1.0	1.2	1.7	0.8	1.4	
	CV	15.6	11.5	13.5	19.2	10.5	17.1	
Galaxy	X	20.5±0.5	21.0 ± 0.6	20.8±0.4	20.1±0.4	19.5±0.4	19.8±0.3	
	s	0.8	0.9	0.8	0.6	0.7	0.7	
	CV	3.7	4.2	4.0	3.0	3.6	3.6	
Gauloises	X	26.0±1.0	25.8±0.8	25.9±0.6	21.6 ± 0.6 1.0 4.6	21.3±0.4	21.4±0.4	
Caporal	s	1.6	1.2	1.4		0.7	0.8	
non-filter	CV	6.2	4.8	5.4		3.2	4.0	
Gauloises Disque Bleu filter	X s CV	21.7 ± 1.2 1.9 8.6	22.7 ± 0.6 1.0 4.4	22.2±0.7 1.6 7.0	17.6 ± 0.5 0.8 4.6	17.1 ± 1.0 1.6 9.2	17.4±0.6 1.2 7.2	

•		<u>''Tar</u>	mg/cigt To	<u>est l</u>	"Tar" mg/cigt Test 2			
BRAND		Smoking Machine A	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine</u> B	A & B Combined	
Half & Half	X s CV	22.6 ± 0.7 1.1 5.0	23.3 ±0.6 1.0 4.4	22.9±0.5 1.1 4.8	21.5 ± 0.4 0.7 3.2	21.0 ± 0.5 0.8 4.0	21.3±0.4 0.8 3.7	
Helmar	X S CV	24.3±0.6 0.9 3.6	24.8±1.0 1.5 6.1	24.6±0.6 1.2 5.0	20.2 ± 0.8 1.3 6.2	20.9±0.5 0.8 3.8	20.6±0.5 1.1 5.2	
Herbert Tareyton	X s CV	29.8±0.9 1.4 4.7	30.7±0.9 1.4 4.4	30.2 ±0.6 1.4 4.8	25.5 ± 0.6 0.9 3.7	25.1±0.5 0.7 2.9	25.3±0.4 0.8 3.3	
Holiday king size, non-filter	X S CV	35.1±1.0 1.6 4.5	33.8 ± 0.7 1.0 3.1	34.4±0.7 1.5 4.3	29.8±0.6 0.9 3.0	28.5±0.5 0.8 3.0	29.2±0.5 1.1 3.7	
Holiday king size, filter	X s CV	26.0±0.9 1.5 5.7	25.3 ±1.1 1.8 7.1	25.7±0.7 1.6 6.4	22.2 ±0.4 0.6 2.9	22.2 ± 0.4 0.7 3.0	22.2±0.3 0.6 2.9	
Home Run	X S CV	20.2 ±0.6 1.0 4.9	20.8±0.7 1.1 5.1	$ \begin{array}{c} 20.5 \pm 0.5 \\ 1.0 \\ 5.2 \end{array} $	17.0±0.4 0.7 4.2	16.1±0.6 1.0 5.9	16.5 ± 0.4 1.0 5.8	
Kent reg. size, filter	X CV	11.1 ± 1.0 1.6 14.8	11.5 ± 1.0 1.5 13.1	11.3±0.7 1.5 13.7	10.0±0.7 1.1 10.9	9.5±0.7 1.2 12.4	9.8±0.5 1.1 11.6	
Kent king size, filter hard pack	X S CV	16.9±0.4 0.7 4.0	17.3±0.5 0.8 4,5	17.1±0.3 0.7 4.3	15.6±1.0 1.5 9.6	15.4±0.4 0.7 4.6	15.5 ±0.5 1.2 7.4	

•		"Tar"	mg/cigt Te	st 1	"Tar	'mg/cigt Te	est 2
BRAND		Smoking Machine A	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine</u> B	A & B Combined
Kent king size, filter soft pack	X S CV	16.1±0.5 0.8 5.1	16.6 ± 0.4 0.7 4.2	16.3±0.3 0.8 4.7	15.4± 0.3 0.5 3.3	15.5 ± 0.6 1.0 6.2	15.4±0.3 0.7 4.8
Kent 100 mm. filter	CV s	19.3±1.0 1.5 7.8	19.7±0.5 0.8 4.2	$\begin{array}{c} 19.5 \pm 0.5 \\ 1.2 \\ 6.2 \end{array}$	19.3±0.3 0.5 2.7	18.9±0.5 0.8 4.2	19.1±0.3 0.7 3.6
King Sano king size, filter	X s CV	6.4±0.7 1.1 17.7	6.6 ± 0.6 1.0 15.4	6.5±0.5 1.0 16.2	7.0 ± 1.2 1.8 25.9	5.9±0.8 1.2 20.1	6.5±0.7 1.6 24.8
King Sano king size, filter menthol	X s CV	7.7 ± 1.0 1.6 20.6	7.2 ± 0.8 1.3 18.3	7.4±0.6 1.4 19.4	7.6 ± 1.4 2.2 28.6	7.0 ± 1.5 2.3 33.2	7.3±1.0 2.2 30.2
Kool reg. size, non-filter, menthol	X S CV	21.5±0.7 1.1 5.1	22.1 ± 0.4 0.7 3.0	21.8± 0.4 0.9 4.3	18.2 ± 0.4 0.6 3.2	18.0±0.8 1.3 7.3	18.0±0.4 1.0 5.5
Kool king size, filter menthol	X S CV	18.3±0.9 1.4 7.5	19.0±0.4 0.6 3.1	18.6±0.5 1.1 5.9	18.1±0.6 1.0 5.6	17.5 ± 0.4 0.7 4.1	17.8±0.4 0.9 5.1
Kool 100 mm. filter, menthol	X S CV	20.2± 0.6 1.0 4.9	20.9 ± 0.5 0.8 3.8	20.6±0.4 1.0 4.6	20.8±0.4 0.6 2.7	20.8±0.4 0.7 3.3	20.8 <u>+</u> 0.3 0.6 2.9

		. "Tar"	"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined		
L & M	X	16.5 ± 0.8	16.7±0.8	16.6±0.6	13.8±0.4	13.5± 0.3	13.7±0.2		
reg. size,	s	1.3	1.3	1.3	0.7	0.4	0.6		
filter	CV	8.1	7.6	7.6	5.0	3.1	4.2		
L & M king size, filter hard pack	X s CV	17.8 ± 0.4 0.7 3.8	17.8±0.6 1.0 5.6	17.8±0.4 0.8 4.6	16.5 ± 0.5 0.8 4.7	16.0±0.5 0.8 4.9	16.2±0.4 0.8 4.9		
L & M king size, filter soft pack	x s CV	$19.7 \pm 0.7 \\ 1.1 \\ 5.4$	19.4±0.7 1.2 6.0	19.5±0.5 1.1 5.6	18.0±0.4 0.6 3.4	18.0±0.6 0.9 5.0	18.0±0.3 0.7 4.2		
L & M	X	20.3±0.9	19.8±0.5	20.0±0.5	19.8±0.8	19.1±0.7	19.4±0.5		
100 mm.	s	1.4	0.8	1.2	1.3	1.1	1.2		
filter	CV	7.1	4.1	5.8	6.4	5.7	6.2		
L & M 100 mm. filter menthol	X s CV	19.4 ± 0.4 0.6 3.2	19.4±0.4 0.6 2.9	19.4±0.2 0.6 3.0	19.3± 0.6 0.9 4.5	19.8±0.6 1.0 5.1	19.6±0.4 1.0 4.9		
Lark	x	18.3±0.6	18.2±0.6	18.3±0.4	16.7±0.8	16.5 ± 0.5 0.8 4.6	16.6±0.4		
king size,	s	0.9	0.9	0.9	1.2		1.0		
filter	CV	4.9	4.8	4.7	7.2		5.9		
Lark	X	18.8±0.5	18.4±0.3	18.6 ±0.3	18.5±0.6	18.4±0.4	18.5 ±0.4		
100 mm.	s	0.8	0.4	0.7	0.9	0.7	0.8		
filter	CV	4.4	2.5	3.6	5.1	3.8	4.4		
Life	X	10.9±0.4	10.8±0.5	10.8±0.3	10.7±0.3	10.9±0.5	10.8±0.3		
	s	0.7	0.8	0.7	0.5	0.8	0.7		
	CV	6.1	7.2	6.5	4.3	7.7	6.1		

		''Tar'	mg/cigt Te	est 1	"Tar	mg/cigt T	st 2	
BRAND		Smoking Machine A	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	
Lucky Strike reg. size, non-filter	X s CV	28.5± 0.6 1.0 3.6	28.0±1.0 1.6 5.7	28.2±0.6 1.3 4.7	23.3±0.6 0.9 4.0	23.0±0.5 0.8 3.4	23.1±0.4 0.8 3.7	
Lucky Strike king size, filter	X s CV	20.9±0.7 1.1 5.2	21.0± 0.4 0.6 2.7	20.9±0.4 0.8 4.0	19.5±0.4 0.7 3.6	19.0±0.5 0.8 4.2	19.3±0.3 0.8 4.0	
Lucky Strike king size, filter menthol	X s CV	19.7 ± 0.7 1.1 5.8	19.7± 0.5 0.8 4.3	19.7±0.4 1.0 5.0	18.2±0.5 0.8 4.6	18.5 ± 0.4 0.7 3.9	18.3±0.3 0.8 4.2	
Lucky Strike 100 mm. filter	X s CV	20.5±0.6 1.0 4.6	20.5 ± 0.6 0.9 4.6	20.5±0.4 0.9 4.5	20.9 ± 0.5 0.8 4.0	20.8±0.5 0.8 3.7	$ \begin{array}{c} 20.8 \pm 0.3 \\ 0.8 \\ 3.7 \end{array} $	
Lucky Strike 100 mm. filter menthol	X s CV	18.2 ± 0.5 0.8 4.3	18.5 ± 0.7 1.1 6.1	18.3 ± 0.4 1.0 5.2	18.9 ± 0.4 0.6 3.4	18.2 ± 0.2 0.3 1.5	18.6±0.2 0.6 3.1	
Mapleton reg. size, non-filter	X S CV	23.9 ± 0.4 0.7 3.0	24.0±0.6 1.0 4.2	23.9±0.4 0.8 3.6	19.4 ± 0.7 1.0 5.4	18.6 ±0.5 0.8 4.1	$ \begin{array}{c} 19.0 \pm 0.4 \\ 1.0 \\ 5.1 \end{array} $	
Mapleton king size, filter	x s cv	21.4 ± 1.0 1.6 7.3	21.1±0.4 0.6 2.8	21.3 ± 0.5 1.2 5.5	20.2±0.4 0.6 3.0	20.0±0.6 0.9 4.4	20.1±0.3 0.7 3.7	
Marlboro king size, filter, hard pack	X S CV	20.6 ± 0.3 0.5 2.6	19.8 ± 1.3 2.0 10.2	20.2 ±0.7 1.5 7.5	18.5±0.3 0.5 2.8	17.9±0.3 0.4 2.5	18.2±0.2 0.6 3.2	

•		<u>''Tar</u>	" mg/cigt Te	est 1	"Tar" mg/cigt Test 2			
BRAND		Smoking Machine A	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	
Marlboro king size, filter soft pack	X s CV	21.7±0.6 1.0 4.5	21.0±0.5 0.9 4.1	21.4±0.4 1.0 4.5	20.2±0.6 0.9 4.4	19.6±0.4 0.6 2.9	19.9±0.3 0.8 3.9	
Marlboro 100 mm. filter soft pack	X S CV	22.4 ± 0.4 0.6 2.9	21.9 ± 0.3 0.5 2.2	22.2 ±0.3 0.6 2.8	22.2 ± 0.4 0.7 3.3	21.4±0.9 1.5 6.9	21.8 ±0.5 1.2 5.5	
Marlboro 95 mm. filter hard pack	X S CV	21.6 ±0.6 1.0 4.8	21.2 ± 0.6 1.0 4.9	21.4±0.4 1.0 4.8	21.6±0.6 1.0 4.8	21.0±0.5 0.8 4.0	21.3±0.4 1.0 4.6	
Marlboro king size, filter , menthol	x s cv	20.3 ± 0.4 0.6 3.0	20.8 ± 0.2 0.4 1.8	20.6 ± 0.2 0.5 2.7	19.6±0.6 0.9 4.6	19.6±0.4 0.6 3.3	19.6±0.3 0.8 3.9	
Marvels reg. size, non-filter	X CV	19.1±0.5 0.8 4.0	19.2±0.6 0.9 4.5	19.2 ± 0.4 0.8 4.2	15.4±0.5 0.8 5.5	15.0±0.4 0.6 3.7	15.2±0.3 0.7 4.7	
Marvels reg. size, filter	X S CV	2.6±0.4 0.6 24.7	3.1±0.3 0.4 14.1	2.8 ±0.3 0.6 20.6	2.8±0.3 0.4 15.9	2.9±0.3 0.5 18.7	2.8±0.2 0.5 17.1	
Marvels king size, non-filter	X S CV	22.8±0.6 1.0 4.4	23.7 ± 1.0 1.7 7.0	23.2 ± 0.6 1.4 6.1	20.0±0.8 1.2 6.2	19.8 ±0.3 0.5 2.5	19.9±0.4 0.9 4.7	
Marvels king size, filter	x cv	5.9±0.4 0.7 12.0	6.8±0.3 0.5 7.3	6.3 ≠0.3 0.8 11.9	5.9±0.4 0.6 11.0	5.7±0.3 0.5 8.4	5.8±0.2 0.6 9.8	

_		"Tar"	mg/cigt Te	<u>st 1</u>	"Tar" mg/cigt Test 2			
BRAND		Smoking Machine A	Smoking Machine B	A & B Combined	Smoking Machine A	Smoking Machine B	A & B Combined	
Marvels	X	4.4±0.2	4.3 ± 0.4	4.3±0.2	4.3 ±0.3	4.2 ± 0.7	4.3±0.4	
king size,	s	0.3	0.7	0.5	0.4	1.0	0.8	
filter menthol	CV	7.4	15.4	11.7	10.2	24.7	18.3	
Montclair	X	9.6 ±0.7	9.8±0.6	9.7±0.4	9.3±0.5	8.7 ± 0.6	9.0±0.4	
	s	1.0	0.9	1.0	0.8	0.9	0.9	
	CV	10.9	9.2	9.9	8.3	10.1	9.6	
Newport king size, filter menthol hard pack	X s . CV	20.5±1.0 1.6 8.0	21.6±0.4 0.6 2.7	21.1±0.6 1.3 6.3	17.6±0.6 1.0 5.4	17.5 ± 0.5 0.8 4.4	17.5±0.4 0.8 4.8	
Newport king size, filter menthol soft pack	X s . CV	20.4± 0.7 1.1 5.6	20.9±0.6 0.9 4.4	$\begin{array}{c} 20.6 \pm 0.5 \\ 1.0 \\ 5.1 \end{array}$	19.0±0.6 1.0 5.1	19.2.±0.4 0.7 3.5	19.1±0.4 0.8 4.2	
Newport 100 mm. filter menthol	X s . CV	20.2±0.9 1.4 6.8	20.4 ± 0.3 0.5 2.3	20.4 ± 0.4 1.0 5.0	20.7 ± 0.6 0.9 4.5	20.2 ± 0.5 0.8 3.8	20.4 ±0.4 0.9 4.2	
Oasis	X	18.9±0.6	19.4±0.6	19.2 ±0.4	17.8 ± 0.5	17.8 ±0.5	17.8±0.3	
	S	1.0	1.0	1.0	0.7	0.8	0.8	
	CV	5.4	5.2	5.4	4.1	4.6	4.3	
Old Gold	X	23.9±0.8	23.3±0.8	23.6 ± 0.6	18.8 ± 0.4 0.6 3.2	18.7±0.4	18.8±0.3	
reg. size	S	1.3	1.4	1.3		0.6	0.6	
non-filter	CV	5.3	5.8	5.6		3.5	3.3	
Old Gold	X	29.2±0.8	29.5±0.8	29.4 ±0.6	24.8 ± 0.7 1.1 4.5	24.4±1.0	24.6±0.6	
king size,	s	1.2	1.3	1.2		1.6	1.4	
non-filter	CV	4.3	4.4	4.2		6.8	5.6	

•		<u>"Tar"</u>	mg/cigt Tes	<u>t 1</u>	"Tar"	t 2	
BRAND		Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B . Combined
Old Gold king size, filter	X s CV	19.5±0.9 1.4 7.2	20.5±1.0 1.5 7.4	20.0±0.7 1.5 7.6	18.8±0.6 1.0 5.5	18.4 ± 0.4 0.6 3.5	18.6±0.4 0.9 4.6
Omar	X S CV	29.7±1.0 1.6 5.2	29.3±0.6 1.0 3.3	29.5±0.6 1.3 4.3	23.7±0.7 1.1 4.7	23.3±0.6 1.0 4.3	23.5±0.5 1.0 4.5
Pall Mall king size, non-filter soft pack	X s CV	28.1±0.9 1.4 4.9	28.4±0.8 1.2 4.2	28.2±0.6 1.3 4.5	23.7±0.6 0.9 3.7	24.1±0.6 0.9 3.9	23.9±0.4 0.9 3.8
Pall Mall 95 mm. filter hard pack	x s CV	18.4±0.4 0.6 3.5	18.1±0.5 0.8 4.7	18.2±0.3 0.7 4.1	18.7 ± 0.5 0.8 4.7	18.3±0.2 0.4 2.1	18.5±0.3 0.7 3.6
Pall Mall 95 mm. filter menthol hard pack	x s CV	17.9±0.6 1.0 5.5	18.1±0.4 0.7 3.9	18.0±0.4 0.8 4.7	18.1±0.5 0.7 4.1	18.3±0.3 0.5 2.8	18.2±0.3 0.6 3.4
Pall Mall 100 mm. filter soft pack	x s CV	20.3±0.4 0.6 2.8	19.7±0.3 0.5 2.7	20.0 ±0.3 0.6 3.1	20.0±0.4 0.6 3.0	19.2±0.5 0.8 4.3	19.6±0.4 0.8 4.2
Pall Mall 100 mm. filter menthol soft pack	X S CV	19.1±0.4 0.7 3.4	18.9±0.4 0.7 3.8	$\begin{array}{c} 19.0 \pm 0.3 \\ 0.7 \\ 3.6 \end{array}$	19.2±0.5 0.8 4.3	18.9±0.4 0.6 3.0	19.0±0.3 0.7 3.7

•		"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	
Parliament king size, filter hard pack	X s CV	15.6 ± 0.5 .0.8 5.3	15.6±0.4 0.7 4.3	15.6 立 0.3 0.7 4.7	15.3±0.5 0.9 5.6	14.7±0.5 0.8 5.2	15.0 ± 0.4 0.8 5.6	
Parliament king size, filter soft pack	x s CV	16.1±0.5 0.8 5.2	16.5 ±0.4 0.6 3.9	16.3±0.3 0.7 4.6	16.5±0.3 0.5 3.3	15.8±0.3 0.5 3.2	16.1±0.3 0.6 3.7	
Philip Morris reg. size, non-filter	X S CV	24.5±0.7 1.1 4.4	25.3±0.5 0.8 3.0	24.9 ± 0.4 1.0 4.0	20.4±0.4 0.6 2.9	20.7±0.6 0.9 4.5	20.6±0.3 0.8 3.7	
Philip Morris Commander	X s CV	31.9±0.8 1.3 4.2	30.6 ± 1.1 1.7 5.6	31.2±0.7 1.6 5.3	26.9 ±1.0 1.6 5.9	26.7±0.5 0.8 3.0	26.8 ±0.6 1.2 4.6	
Philip Morris king size, filter	X s CV	18.8±0.5 0.9 4.6	18.8±0.6 0.8 4.6	18.8±0.4 0.8 4.5	18.5±0.3 0.5 2.9	18.4±0.5 0.8 4.7	18.4±0.3 0.7 3.8	
Philip Morris king size, filter menthol	X s L CV	19.9 ±0.6 1.0 5.1	20.0±0.8 1.2 6.0	19.9 ± 0.5 1.1 5.4	19.1±0.5 0.8 4.4	19.2±0.5 0.8 4.0	19.1±0.3 0.8 4.1	
Picayune	X S CV	20.3±0.5 0.8 3.7	20.5±0.3 0.5 2.4	20.4±0.3 0.6 3.1	16.8±0.6 1.0 6.2	16.3 ±0.4 0.6 3.8	16.6±0.4 0.9 5.2	

		"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking Machine	Smoking Machine	A & B Combined	Smoking Machine	Smoking Machine	A & B Combined	
Piedmont	X	24.7±0.8	25.1±0.6	24.9±0.5	21.1±0.4	20.5±1.0	20.8 ±0.5	
	S	1.3	0.9	1.1	0.6	1.5	1.2	
	CV	5.1	3.7	4.4	3.1	7.4	5.7	
Players	X S CV	32.1 ± 0.8 1.2 3.8	31.5±0.5 0.8 2.6	31.8±0.5 1.1 3.3	26.2±0.5 0.7 2.8	25.4±0.5 0.8 3.1	25.8±0.4 0.8 3.3	
Raleigh	X	30.0 ± 0.8 1.3 4.2	29.8±1.2	29.8±0.7	25.5±0.4	25.2 [±] 0.6	25.4 ±0.4	
king size,	s		1.8	1.5	0.6	1.0	0.8	
non-filter	CV		6.1	5.2	2.4	3.9	3.3	
Raleigh	X	18.4±0.6	19.1±0.7	18.7±0.5	17.6±0.4	17.8±0.5	17.7±0.3	
king size,	S	0.9	1.1	1.0	0.6	0.8	0.6	
filter	CV	5.0	5.6	5.5	3.2	4.3	3.7	
Salem	X	19.9±0.5	20.5 ± 0.4	20.2 ±0.4	18.5±0.5	18.6±0.3	18.6 ±0.3	
king size,	s	0.9	0.7	0.8	0.8	0.4	0.6	
filter menthol	CV	4.3	3.2	4.0	4.3	2.2	3.3	
Salem	X	21.3± 0.4	21.5±0.8	21.4±0.4	21.1±0.4	21.2±0.5	21.1±0.3	
100 mm.	s	0.6	1.3	1.0	0.6	0.7	0.7	
filter menthol	. CV	2.6	6.0	4.6	2.7	3.5	3.3	
Sano	X	15.5±1.0	15.3±1.0	15.4±0.7	12.5±0.6	12.3±0.7	12.4±0.4	
reg. size,	S	1.6	1.6	1.6	0.9	1.1	1.0	
non-filter	CV	10.8	10.7	10.5	7.0	8.9	7.8	
Sano	X	3.1±0.4	3.6±0.3	3.4 ±0.3	3.5±0.2	3.1±0.4	3.3 ± 0.2 0.5 16.0	
reg. size,	S	0.6	0.5	0.6	0.4	0.6		
filter	CV	21.0	14.0	18.5	10.9	18.7		

	•		'Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
	BRAND		Smoking <u>Machine A</u>	Smoking Machine B	A & B Combined	Smoking Machine A	Smoking Machine B	A & B Combined	
	Silva Thins	X	15.1±0.5	15.0±0.4	15.0±0.3	14.8±0.5	14.5±0.3	14.7±0.3	
	100 mm.	s	0.8	0.6	0.7	0.8	0.5	0.6	
	filter	CV	5.1	4.3	4.6	5.2	3.5	4.4	
	Silva Thins	X	14.6±0.5	14.7±0.5	14.6±0.3	15.0 ±0.3	15.1±0.5	15.1 ±0.3	
	100 mm.	s	0.7	0.8	0.7	0.5	0.8	0.7	
	filter menthol	CV	5.0	5.4	5.1	3.5	5.7	4.6	
	Spring	X s CV	20.5 ± 0.8 1.3 6.5	20.5±0.7 1.1 5.6	20.5±0.5 1.2 5.9	21.0 ±0.7 1.1 5.4	20.9±0.6 0.9 4.2	21.0±0.4 1.0 4.7	
	Stratford	X	25.3±0.8	25.4 ±0.9	25.3±0.6	21.3 ±0.6	21.0±0.3	21.2±0.3	
	king size,	s	1.2	1.4	1.2	0.9	0.5	0.7	
	non-filter	CV	4.8	5.4	5.0	4.1	2.6	3.4	
	Stratford	X	18.7 <u>十</u> 0.5	18.6±0.5	18.7±0.4	17.5 ±0.3	17.2±0.2	17.4±0.2	
	king size,	s	0.8	0.8	0.8	0.5	0.4	0.4	
	filter	CV	4.4	4.5	4.4	2.7	2.2	2.5	
	Stratford	X	19.7±0.6	20.1±0.4	19.9±0.4	18.5 ±0.5	18.9±0.4	18.7±0.3	
	king size,	s	0.9	0.7	0.8	0.8	0.6	0.7	
	filter menthol	CV	4.5	3.3	4.0	4.4	3.0	3.8	
	Sweet Caporal	X s CV	21.6±0.9 1.4 6.4	22.1±0.5 0.8 3.8	21.8 ± 0.5 1.2 5.3	20.9±0.6 1.0 4.8	19.9±0.5 0.8 4.0	20.4 <u>十</u> 0.5 1.0 5.1	
	Tareyton	X	16.8±0.7	16.5±0.4	16.6±0.4	15.8 ± 0.2	15.6±0.5	15.7±0.3	
	king size,	s	1.0	0.6	0.8	0.3	0.9	0.6	
	filter	CV	6.3	3.5	5.1	1.6	5.5	4.0	
1	Tareyton	X	17.2±0.6	16.9±0.3	17.1±0.3	17.0±0.4	16.8±0.5	17.0±0.3	
	100 mm.	s	0.9	0.5	0.7	0.6	0.8	0.7	
	filter	CV	5.4	3.0	4.4	3.7	5.0	4.3	
					,				

		"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	
Tempo	X	13.5±0.8	13.5±0.9	13.5±0.6	13.4±0.6	12.8±0.3	13.1±0.4	
	s	1.2	1.4	1.3	1.0	0.5	0.8	
	CV	9.2	10.1	9.4	7.5	3.9	6.3	
True (Blue)	X	12.3±0.6	13.0 ±0.5	12.6 ±0.4	12.8±0.3	12.6 ± 0.3 0.5 3.9	12.7±0.2	
king size,	s	1.0	0.8	0.9	0.5		0.5	
filter	CV	8.0	6.2	7.5	4.0		3.9	
True (Green)	X	12.5 ± 1.2 2.0 15.8	12.2±0.4	12.4±0.6	12.9±0.6	12.5±0.4	12.7 ±0.3	
king size,	s		0.7	1.4	0.9	0.6	0.8	
filter menthol	CV		.5.7	11.7	7.0	4.9	6.1	
Tryon	X	10.8±0.7	11.2 ±0.5	11.0 ±0.4	10.2±0.3	10.3±0.3	10.2±0.2	
king size,	s	1.0	0.7	0.9	0.5	0.4	0.5	
filter	CV	9.6	6.6	8.1	5.0	4.5	4.6	
Tryon	X	11.2±1.0	11.0±0.6	11.1±0.6	10.1±0.4	10.1±0.3	10.1±0.3	
king size,	s	1.5	1.0	1.2	0.6	0.6	0.6	
filter menthol	C V	13.6	8.9	11.2	6.5	5.5	5.9	
Viceroy	X	18.9±0.5	19.0±0.5	$ \begin{array}{c} 19.0 \pm 0.3 \\ 0.8 \\ 4.1 \end{array} $	17.4±0.5	17.8±0.5	17.6±0.3	
king size,	S	0.8	0.8		0.8	0.7	0.8	
filter	CV	4.2	4.0		4.6	4.2	4.4	
Viceroy	X	19.6±0.7	19.9±0.5	19.7 ±0.4	20.0±0.8	20.1±0.4	20.0±0.4	
100 mm.	s	1.2	0.8	1.0	1.2	0.7	0.9	
filter	CV	5.9	4.1	5.0	6.0	3.3	4.7	
Virginia Slims	X	17.9 ±0.7	18.7±0.7	18.3±0.5	19.0±0.4	18.2±0.5	18.6±0.4	
100 mm.	s	1.2	1.1	1.2	0.7	0.8	0.8	
filter	CV	6.5	6.0	6.5	3.5	4.4	4.3	

		"Tar" mg/cigt Test 1			"Tar" mg/cigt Test 2			
BRAND		Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	
Virginia Slims 100 mm. menthol	x s CV	18.3±0.5 0.8 4.1	18.7±0.5 0.9 4.6	18.5±0.4 0.8 4.5	19.2±0.3 0.5 2.8	18.2±0.3 0.4 2.4	18.7±0.3 0.7 3.7	
Vogue Assorted Colors	X s CV	11.6±0.5 0.8 7.2	11.6±0.6 1.0 8.3	11.6±0.4 0.9 7.6	8.8±0.4 0.6 6.4	9.1±0.3 0.5 5.4	8.9±0.2 0.5 6.0	
Vogue black	X s CV	18.4±1.1 1.8 9.7	19.4 ± 1.6 2.6 13.2	18.9±1.0 2.2 11.6	15.4±0.7 1.1 7.0	15.3±0.8 1.2 7.8	15.3±0.5 1.1 7.2	
Wings	X s CV	19.7±0.5 0.8 4.2	19.4±0.6 0.9 4.5	19.6 ±0.4 0.8 4.4	18.2±0.4 0.6 3.4	18.3±0.6 0.9 5.2	18.2 ±0.3 0.8 4.3	
Winston king size, filter hard pack	x s CV	21.1±0.4 0.7 3.4	21.0±0.6 1.0 4.6	21.1±0.4 0.8 3.9	17.4±0.4 0.7 4.0	17.5±0.2 0.4 2.3	17.5±0.2 0.6 3.2	
Winston king size, filter soft pack	X s CV	20.5±0.4 0.7 3.4	20.9±0.7 1.1 5.2	20.7±0.4 0.9 4.4	19.1±0.3 0.5 2.8	19.1±0.6 0.9 4.6	19.1±0.3 . 0.7 3.7	
Winston 100 mm. filter, soft pack	x s CV	21.9 ± 0.4 0.7 3.3	21.4±0.5 0.8 3.6	21.6±0.3 0.8 3.5	21.4±0.3 0.5 2.2	21.3±0.4 0.6 2.9	21.4±0.2 0.5 2.5	

U	
Ž	
$\overline{\omega}$	
9	
\simeq	
8	
Ğ	
N	
\sim	
$\frac{\omega}{\Delta}$	

•		Tar	mg/cigt To	est <u>l</u>	"Tar" mg/cigt Test 2		
BRAND		Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined	Smoking <u>Machine A</u>	Smoking <u>Machine B</u>	A & B Combined
Winston 100 mm. filter menthol	X S CV	21.5±0.6 0.9 4.3	21.6±0.4 0.7 3.0	21.6±0.3 0.8 3.6	21.3 ±0.6 0.9 4.2	20.9±0.6 0.9 4.4	21.1±0.4 0.9 4.3
Yukon	X s CV	19.3±0.8 1.3 6.6	19.6±0.6 1.0 5.0	19.4±0.5 1.1 5.7	16.2±0.4 0.7 4.2	15.7±0.9 1.5 9.5	15.9±0.5 1.2 7.3
Monitor Smoking Machine A	X s CV	19.9 <u>†</u> 0.1 1.1 5.4	19.9±0.1 1.0 5.2		20.0±0.1 0.7 3.5	$\begin{array}{c} 19.8 \pm 0.1 \\ 0.7 \\ 3.3 \end{array}$,